

Nephchem Complement C3 (Nephelometry method)



KIT NAME	KIT SIZE	CAT. NO
Nephchem - C3	25 Tests	NCC301025T

INTRODUCTION

Complement C3 (C3) is intended for Invitro quantitative determination of C3 in human serum. Complement component 3 (C3) is the most abundant protein of the complement system. It plays a central role in the complement system and contributes to innate immunity. When activated by classical or alternative pathway, complement destroys foreign agents, release histamine, and release leukocytes from the bone marrow. C3 levels can be used in determining inherited or acquired deficiencies, including active lupus nephritis, severe infections, and inflammation. Increased levels are usually found after trauma, surgery, biliary obstruction, or focal glomerulosclerosis.

METHOD PRINCIPLE

The kit utilizes latex-enhanced immunoturbidimetry to measure the C3 level in human serum by GB NEPHCHEM (Nephelometry method). The Kit utilizes immunoturbidimetry to measure the C3 level in human serum. During the test, C3 in the sample binds with the specific anti-C3 antibody to cause agglutination. The turbidity caused by agglutination is detected optically by GB NEPHCHEM, analyzer. The change in absorbance is proportional to the level of C3 in the sample. The actual concentration is obtained by comparing with a calibration curve with known concentrations.

KIT CONTENTS

Reagent kit - box	
R1 - C3 buffer	1x 4.9 ml
R2 - C3 antibody	1 x 1.6 ml
Test Card	1 no
Accessories kit box	
Cuvettes	25 nos
Big tips	25 nos
small tips	50 nos

Working reagent preparation and stability

Reagent R1 and R2 are ready to use liquid stable at 2-8°C till the expiry date printed on the package.

Concentrations in the test

R1 - Phosphate buffer, Polyethylene glycol, Sodium azide < 0.1%
R2 - anti-C3 antibodies, Tris buffer, sodium azide < 0.1%

Warnings and notes

1. The Kit is for *in vitro* diagnostic use only. Not for use in humans or animals.
2. The instructions must be followed to obtain accurate results.
3. Do not use the reagents beyond the expiration date.
4. Treat all specimens as infectious. Proper handling and disposal procedures of specimens and test materials should be strictly followed.
5. Reagents contain less than 0.1% sodium azide as preservative; avoid contact with skin and eyes, flush with copious amounts of water when disposing.

SPECIMEN

Follow standard laboratory procedures to collect serum samples. It is recommended to perform test immediately after sample collection. If the test cannot be done immediately, store sample at 2-4°C for up to 3 days or at -20°C for up to 1 months. Avoid repeated freezing and thawing.

PROCEDURE

It is very important for antigen-antibody reaction needs the pre-warm of both reagents and samples. Along with GB NEPHCHEM equipment, dry bath incubator will be provided, in that dedicated R1, R2 and sample positions were available. Please use the respective positions for desired pre-warm temperature of 37°C

- Step 1: Insert Test Card to Card reader slot and display will show promptly add R1 + S (sample)
- Step 2: Pipette out 180 µl of R1 into dedicated cuvette and add 5 µl of sample (serum) and place the cuvette in the reading chamber
- Step 3: After the incubation, the display will show promptly add R2
- Step 4: Pipette out 60 µl of R2 using sensor pipette connected with machine into the cuvette
- Step 5: Once the reaction time got over, the result will show in the display and (if external printer connected then it will get print out)

REFERENCE VALUES

80 to 185 mg/dL

It is recommended for each laboratory to establish its own reference ranges for local population.

QUALITY CONTROL

To ensure adequate quality control, each kit can be cross checked with commercially available third party Immunological quality control or use recommended GB Immunology Quality control.

PERFORMANCE CHARACTERISTICS

- **Linearity:** 0 to 335.0 mg/dL
- **Precision:** within Run CV ≤ 6 %
- **Specificity / Interferences**
No interference detected for bilirubin upto 60 mg/dL and hemoglobin 10 g/L, triglycerides 1000 mg/dL

WASTE MANAGEMENT

Please refer to local legal requirements.

LITERATURE

1. Burtis C., Ashwood, E.R.(ed). Tietz Textbook of Clinical Chemistry, 3rd ed. Philadelphia, PA; WB Saunders Co; 503-4; 1999.
2. De Bruijn MH, Fey GH (February 1985). *oc. Natl. Acad. Sci. U.S.A.* 82 (3): 708-12.
3. Lachmann P (December 1975). *Med. Genet.* 12 (4): 372-7.
4. Matsuyama W, Nakagawa M, Takashima H, Muranaga F, Sano Y, Osame M (December 2001). "Molecular analysis of hereditary deficiency of the third component of complement (C3) in two sisters". *Intern. Med.* 40 (12): 1254-8.
5. Sahu A, Lambris JD (April 2001). "Structure and biology of complement protein C3, a connecting link between innate and



Genuine Biosystem Private Limited
Plot No.97 & 98, kattabomman street,
Parvathy Nagar Extension,
Old Perungalathur, Chennai - 600063, India.
Ph: +91-44-48681845
Email: genuinebiosystem@gmail.com
website: www.genuinebiosystem.com