



Genuine Biosystem Pvt Ltd

Genuine Service for better tomorrow

PHOTIC 10

DIGITAL



PHOTOELECTRIC COLORIMETER



2½ —

DIGIT LED DISPLAY



1.5 —

KG WEIGHT



8 —

OPTICAL GLASS FILTERS

PHOTIC 10

PHOTOELECTRIC COLORIMETER

DIGITAL

TECHNICAL SPECIFICATIONS

Filters	8 (All High Quality Optical Glass Filters)
Filters Range	400,480,510,520,540,580,620,700 (nm)
Min. Volume	1ml.
Display	2½ digit LED
Range	400-700nm
Output	OD: 0 to 1.99
Resolution	OD: 0.01
Accuracy	0.5% FS ± 1 Digit
Detector	Hermetically Sealed Imported Photodiode
Light Source	Imported White LED
Power	230V ± 10% AC, 50 Hz
Dimensions	L-33cm. x W-30cm. x H-16.5cm. (approx.)
Weight	1.5 kg. Approx.
Accessories	(i) Test Tube-5 (ii) Test Tube Stand (iii) Dust Cover-1 (iv) Spare Light Source-2 (v) Instruction Manual-1

Technical Specifications Subject To Change Without Prior Notice.



Genuine Biosystem Pvt Ltd

Genuine Service for better tomorrow



PRODUCTION SITE :

Plot no. 97 & 98, Kattabomman street, Parvathy Nagar Extension, Old Perungalathur, Chennai-600063, India.

SALES OFFICE :

No.8/51 , 2nd floor , Kalaignar Nedunchalai, Peerkankaranai , Chennai 600063, India.

📧 sales@gb-group.co.in | marketing@gb-group.co.in

📞 +91 – 9790943607 | 9003039750

📞 +91 – 44 – 48681845