



## **MICROPLATE ELISA WASHER** **ENSURE COMPLETE WASHING**

- Large LCD display, menu designed for convenient use
- Whole plate washing or single strip washing
- 12-way and 8-way manifold included
- Low residual volume by two pipettes
- Complete bottom washing
- Automatic monitoring of vacuum and pressure, automatic rinse cycle
- Flat, V-Bottom or U-Bottom plate and strips washing
- Large memory to store up to 48 user programmed wash protocols
- Emergency stop available and liquid warning



# GB LISA W

## A RAPID MICROPLATE ELISA WASHER



2 Modes of operating strip and plate



2 Point as earlier



100 User defined wash programs

### TECHNICAL SPECIFICATION

#### INPUT

8 Key membrane keyboard

#### WASH MODE

Strip mode and plate mode

#### WASHING TIME

0~99 times adjustable

#### RESIDUAL VOLUME

≤ 1μL/Well

#### STORAGE CAPACITY

More than 100 user defined wash board procedures

#### NET WEIGHT

8kg

#### WASH STRIPS

1~12 adjustable

#### DIPPING TIME

0~3600s adjustable

#### LIQUID VOLUME

50~3000μL/Well

#### POWER SUPPLY

198~242V, 49~51Hz

#### OPERATING ENVIRONMENT

Temperature 5°C~40°C, Max humidity 80%

#### MANIFOLD

8 pins and 12 pins

#### APPLICABLE WELL

Flat, V-Bottom or U-bottom

#### SHAKING TIME

0~600s adjustable

#### SHIPPING TIME

0.1~9.9s

#### DIMENSIONS

448 (L) x 382 (W) x 163 (H) mm



Genuine Biosystem Pvt Ltd

#### Factory

Plot No. 97 & 98, Kattabomman street, Parvathy Nagar extension, Old Perungalathur, Chennai 63, India .

#### Sales

Plot No.54, Kambar Street, 1st & 2nd Floor, Parvathy Nagar, Old Perungalathur, Chennai 600063, TN (India).



sales@genuinebiosystem.com | export@genuinebiosystem.com



+91 9003039750 | 9677126493 | 9677235540



+91 – 44 - 48681845