



Specific Protein Analyzer

- Nephelometry Principal
- Innovative smart card technology
- Easy Automated Calibration
- Simple procedure No reagent wastage
- Auto Rotation / Auto Mixing /Auto Incubation
- 650 nm LED Light Source (Long life)
- Small Pack Size make its convenience
- User Friendly Interface



GB NEPHCHEM

A SPECIFIC PROTEIN ANALYZER

EASY AUTOMATED CALIBRATION



1 SINGLE CHANNEL
LED



100 RESULTS
MEMORY

TECHNICAL SPECIFICATIONS

PRINCIPLE

Nephelometry

LIGHT SOURCE

LASER LED 650nm

DISPLAY

2Line LED Back Light Display

SAMPLE TYPE

Plasma, Serum, Whole Blood, Urine

INCUBATION FOR SAMPLES

16 Position, Separate Dry Bath

INCUBATION FOR REAGENTS

2 Position, R1 and R2

TESTING POSITION

Single Channel

MIXING

Automatic By Motor

PRINTER

External Printer Via Rs232 (Option)

ASSAY TEMPERATURE

37.0 ± 1.0°C

ENVIRONMENTAL CONDITIONS

Temperature 10-30°C

HUMIDITY

20%-70%

WEIGHT

Less Than 1.5 kg

WIDE RANGE OF PARAMETERS

Product Name	Pack Size
ASO	25 T
CRP	25 T
RF	25 T
hSCRp	25 T
HbA1C	25 T
MICROALBUMIN	25 T
IgA	25 T
IgE	25 T
IgM	25 T
IgG	25 T
APO A1	25 T
APO B	25 T
C3	25 T
C4	25 T
Lp (a)	25 T
CYSTATIN C	25 T
D-DIMER	25 T

SMARTCARD CALIBRATION



Genuine Biosystem Pvt Ltd

Genuine Service for better tomorrow



PRODUCTION SITE :

Plot no. 97 & 98, Kattabomman street, Parvathy Nagar Extension, Old Perungalathur, Chennai-600063, India.

SALES OFFICE :

No.8/51 , 2nd floor , Kalaignar Nedunchalai, Peerkankaranai , Chennai 600063, India.

✉ sales@gb-group.co.in | marketing@gb-group.co.in

☎ +91 – 9790943607 | 9003039750 ☎ +91 – 44 – 48681845