



## Specific Protein Analyzer

- Nephelometry Principal
- Innovative smart card technology
- Easy Automated Calibration
- Simple procedure No reagent wastage
- Auto Rotation / Auto Mixing / Auto Incubation
- 650 nm LED Light Source (Long life )
- Small Pack Size make its convenience
- User Friendly Interface



# GB NEPHCHEM

## A SPECIFIC PROTEIN ANALYZER

EASY AUTOMATED CALIBRATION



**1** SINGLE CHANNEL  
LED



**100** RESULTS  
MEMORY

### TECHNICAL SPECIFICATIONS

#### PRINCIPLE

Nephelometry

#### LIGHT SOURCE

LASER LED 650nm

#### DISPLAY

2Line LED Back Light Display

#### SAMPLE TYPE

Plasma, Serum, Whole Blood, Urine

#### INCUBATION FOR SAMPLES

16 Position, Separate Dry Bath

#### INCUBATION FOR REAGENTS

2 Position, R1 and R2

#### TESTING POSITION

Single Channel

#### MIXING

Automatic By Motor

#### PRINTER

External Printer Via Rs232 (Option)

#### ASSAY TEMPERATURE

37.0 ± 1.0°C

#### ENVIRONMENTAL CONDITIONS

Temperature 10-30°C

#### HUMIDITY

20%-70%

#### WEIGHT

Less Than 1.5 kg

#### WIDE RANGE OF PARAMETERS

Product Name	Pack Size
ASO	25 T
CRP	25 T
RF	25 T
hSCRp	25 T
HbA1C	25 T
MICROALBUMIN	25 T
IgA	25 T
IgE	25 T
IgM	25 T
IgG	25 T
APO A1	25 T
APO B	25 T
C3	25 T
C4	25 T
Lp (a)	25 T
CYSTATIN C	25 T
D-DIMER	25 T

#### SMARTCARD CALIBRATION



**Genuine Biosystem Pvt Ltd**

Genuine Service for better tomorrow



#### PRODUCTION SITE :

Plot no. 97 & 98, Kattabomman street, Parvathy Nagar Extension, Old Perungalathur, Chennai-600063, India.

#### SALES OFFICE :

No.8/51 , 2nd floor , Kalaignar Nedunchalai, Peerkankaranai , Chennai 600063, India.

sales@gb-group.co.in | marketing@gb-group.co.in

+91 - 9790943607 | 9003039750 | +91 - 44 - 48681845